


Sizing Guide


For Use With
Simplus Full Face
Vitera Full Face



How to Use:

1. Find a ruler and a mirror
2. Face the mirror and put the ruler perpendicular to your face.
3. Measure from the brow indent to the indent below the lower lip
4. Consult our size guide below.

Results will vary and this sizing guide will NOT guarantee the most optimal fit with the mask.



Example: Medium

Size Guide:

Small:	3.0 - 3.6 Inches	7.5 - 9.0 cm
Medium:	3.6 - 4.0 Inches	9.0 - 10.0 cm
Large:	4.0+ Inches	10.0+ cm

Note: We find these measurements work with about 80% of our clients. Head shape, nose shape and nose density is also a factor. In person fitment is always recommended.

