
Sizing Guide

For Use With Wisp DreamWisp






How to Use:

- 1. Find a ruler and a mirror.
- 2. Face the mirror and put the ruler horizontal to your face.
- 3. Measure from outer nostril to outer nostril.
- 4. Consult our size guide below.

Results will vary and this sizing guide will NOT guarantee the most optimal fit with the mask.



Size Guide:

Petite:	0 - 1.1 Inches	0 - 2.8 cm
Small:	1.1 - 1.3 Inches	2.8 - 3.3 cm
Medium:	1.3 - 1.6 Inches	3.3 - 4.0 cm
Large:	1.6 - 1.75 Inches	4.0 - 4.5 + cm
Extra Large:	1.75 + Inches	4.5 + cm

Note: We find these measurments work with about 80% of our clients. Head shape, nose shape and nose density is also a factor. In person fitment is always recommended.



